


# Compatibility Innovation Value

**BIO | One Stage**



# 2016

**Product Catalog**

**BIO | One Stage**

- Straumann Compatible

Valid January 2016



**BlueSkyBio**

## Table of contents

<b>Implants</b>	<b>2,3</b>
<b>Special Cover Screws</b>	<b>3</b>
<b>Solid Abutments and Restorative Kits</b>	<b>4,5</b>
<b>synOcta® Compatible Abutments</b>	<b>5,6</b>
<b>Gold UCLA Abutments</b>	<b>6,7</b>
<b>Angled Abutments</b>	<b>7</b>
<b>Overdenture Abutments</b>	<b>7,8</b>
<b>Intraorally Preparable Abutments</b>	<b>8,9</b>
<b>Waxing Sleeves / Temporary Transfers</b>	<b>9</b>
<b>Abutment level analogs - Metal</b>	<b>9,10</b>
<b>Impression Transfers and Implant Level Analogs</b>	<b>10</b>
<b>Surgical Instruments</b>	<b>11,12</b>
<b>X Cube, Suture</b>	<b>13</b>

<b>Master Bone Expander Kit</b>	<b>13,14</b>
<b>Sinus Lift Kit, Safety Lateral Sinus Kit</b>	<b>14</b>
<b>Terms &amp; Conditions</b>	<b>15</b>



## Implants

Implants blasted with resorbable media and acid washed - Cover screw included.

Compatible with Straumann® ITI® implants, SynOcta® system, solid abutments, instrumentation and prosthetics. Our new narrow platform uses the same prosthetics as regular platform.

For easy insertion, implants must be placed with a J-implant insertion driver and should not be used with other companies' mounts.

### Narrow 3.3

Length	4.1mm	Size	Part #
		Ø3.3 x 10mm, 4.1mm Platform One Stage Implant	ION3310
		Ø3.3 x 12mm, 4.1mm Platform One Stage Implant	ION3312

### Regular 3.3

Length	4.8mm	Size	Part #
		Ø3.3 x 10mm, Regular Platform One Stage Implant	IOR3310
		Ø3.3 x 12mm, Regular Platform One Stage Implant	IOR3312
		Ø3.3 x 14mm, Regular Platform One Stage Implant	IOR3314

### Regular 4.8

Length	4.8mm	Size	Part #
		Ø4.8 x 8mm, Regular Platform One Stage Implant	IOR4808
		Ø4.8 x 10mm, Regular Platform One Stage Implant	IOR4810
		Ø4.8 x 12mm, Regular Platform One Stage Implant	IOR4812
		Ø4.8 x 14mm, Regular Platform One Stage Implant	IOR4814
		Ø4.8 x 16mm, Regular Platform One Stage Implant	IOR4816

### Narrow 4.1

Length	4.1mm	Size	Part #
		Ø4.1 x 10mm, 4.1mm Platform One Stage Implant	ION4110
		Ø4.1 x 12mm, 4.1mm Platform One Stage Implant	ION4112

### Regular 4.1

Length	4.8mm	Size	Part #
		Ø4.1 x 8mm, Regular Platform One Stage Implant	IOR4108
		Ø4.1 x 10mm, Regular Platform One Stage Implant	IOR4110
		Ø4.1 x 12mm, Regular Platform One Stage Implant	IOR4112
		Ø4.1 x 14mm, Regular Platform One Stage Implant	IOR4114
		Ø4.1 x 16mm, Regular Platform One Stage Implant	IOR4116


### Wide 4.8


Length	6.5mm	Size	Part #
		Ø4.8 x 8mm, Wide Platform One Stage Implant	IOW4808
		Ø4.8 x 10mm, Wide Platform One Stage Implant	IOW4810
		Ø4.8 x 12mm, Wide Platform One Stage Implant	IOW4812


All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.


Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Wide 5.6		
	Size	Part #
	Ø5.6 x 8mm, Wide Platform One Stage Implant	ITW5608
	Ø5.6 x 10mm, Wide Platform One Stage Implant	ITW5610

Wide 5.6		
	Size	Part #
	Ø4.8 x 6mm, Regular Platform One Stage Implant	IORS

Wide 5.6		
	Size	Part #
	Ø5.6 x 6mm, Wide Platform One Stage Implant	IOWS

## Special Cover Screws

Internal Cover Screws		
	Size	Part #
	0mm - Internal Cover Screw, One Stage System	CF0010

1mm Cover Screws		
	Size	Part #
	1mm - Narrow Cover Screw, One Stage System	CN0010
	1mm - Regular Cover Screw, One Stage System	CR0010
	1mm - Wide Cover Screw, One Stage System	CW0010

3mm Cover Screws		
	Size	Part #
	3mm - Narrow Cover Screw, One Stage System	CN0030
	3mm - Regular Cover Screw, One Stage System	CR0030
	3mm - Wide Cover Screw, One Stage System	CW0030

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann



## Solid Abutments and Restorative Kits






Master Restorative Kits include a solid abutment, one metal analog (for finishing restoration), two resin analogs (for model), snap impression transfer, engaging and non engaging waxing sleeves.

Note: For single crowns use the clear engaging sleeve for undercasting and the white sleeve is used for temporary cap. For splinted restorations use the white sleeve for the undercasting and the clear one for the temporary cap.

### Solid Prosthetics Regular Platform 4.0mm

	Size	Part #
	Master Restorative Kit - Regular 4.0mm	KRT0040
	Abutment, regular 4mm	AR4000
	Snap Impression Transfer, Reg. 4mm	MSR4000
	Abutment Level Analog 4mm	BAR4000
	Waxing/Temporary Sleeves Regular	WIRK4800

### Solid Prosthetics Regular Platform 5.5mm

	Size	Part #
	Master Restorative Kit - Regular 5.5mm	KRT0055
	Abutment, regular 5.5mm	AR5500
	Snap Impression Transfer, Reg. 5.5mm	MSR5500
	Abutment Level Analog 5.5mm	BAR5500
	Waxing/Temporary Sleeves Regular	WIRK4800

### Solid Prosthetics Regular Platform 7.0mm

	Size	Part #
	Master Restorative Kit - Regular 5.5mm	KRT0070
	Abutment, regular 5.5mm	AR7000
	Snap Impression Transfer, Reg. 5.5mm	MSR7000
	Abutment Level Analog 5.5mm	BAR7000
	Waxing/Temporary Sleeves Regular	WIRK4800

### Solid Prosthetics Narrow Platform 4.0mm

	Size	Part #
	Master Restorative Kit - Narrow 4.0mm	KNT0040
	Abutment, regular 5.5mm	AR4000
	Snap Impression Transfer, Narrow 5.5mm	MSN4000
	Abutment Level Analog 4mm	BAN4000
	Waxing/Temporary Sleeves Narrow	WINK4100

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

## Solid Prosthetics Narrow Platform 5.5mm

	Size	Part #
	Master Restorative Kit - Narrow 5.5mm	KNT0055
	Abutment, regular 4mm	AR5500
	Snap Impression Transfer, Reg. 4mm	MSN5500
	Abutment Level Analog 4mm	BAN5500
	Waxing/Temporary Sleeves Narrow	WINK4100

## Solid Prosthetics Narrow Platform 7.0mm

	Size	Part #
	Master Restorative Kit - Regular 5.5mm	KNT0070
	Abutment, regular 5.5mm	AR7000
	Snap Impression Transfer, Reg. 5.5mm	MSN7000
	Abutment Level Analog 5.5mm	BASN7000
	Waxing/Temporary Sleeves Regular	WINK4800

## Solid Prosthetics Wide Platform 4.0mm

	Size	Part #
	Master Restorative Kit - Wide 4.0mm	KWT0040
	Abutment, wide 4.0mm	AW4000
	Snap Impression Transfer, Wide 4mm	MSW4000
	Abutment Level Analog 4.0mm	BAW4000
	Waxing/Temporary Sleeves Wide	WIWK6500



## Solid Prosthetics Wide Platform 5.5mm

	Size	Part #
	Master Restorative Kit - Wide 5.5mm	KWT0055
	Abutment, wide 5.5mm	AW5500
	Snap Impression Transfer, Wide 5.5mm	MSW5500
	Abutment Level Analog 5.5mm	BAW5500
	Waxing/Temporary Sleeves Wide	WIWK6500

## synOcta® Compatible Abutments

Abutments include a fixation screw. Only Compatible with Tissue Level implants. Compatible with Straumann®-ITI® Lifecore® Stage 1® implants. synOcta® compatible. Must be inserted with .048 hex driver. Regular abutments can be used with the regular and narrow platform implants.

### synOcta® Compatible Abutments


	Size	Part #
	Abutment Oct Regular/Narrow	AOR5500
	Abutment Oct Wide	AOW5500


All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann


synOcta® Compatible Abutments		
	Size	Part #
	Custom Millable Regular/Narrow	AOP5500
	Custom Millable Wide	AOWP5500


## Gold UCLA Abutments

UCLA abutments are gold alloy.  
Only Compatible with Tissue Level implants.

Gold UCLA Prosthetics Regular Platform Engaging		
	Size	Part #
	Restorative Kit Gold UCLA - Regular	KAUCLAORG
	Gold UCLA Abutment, Regular (including screw)	AUCLAORG
	Impression Transfer - Snap Regular	MPOR
	Analog – Implant Level Regular	BOR48
	Abutment Screw, UCLA, One Stage	STT20UC


Gold UCLA Prosthetics Regular Platform Engaging		
	Size	Part #
	Restorative Kit Gold UCLA - Regular	KAUCLAORNG
	Gold UCLA Abutment Non-Engaging, Regular (including screw)	AUCLORNG
	Impression Transfer - Snap Regular	MPOR
	Analog – Implant Level Regular	BOR48
	Abutment Screw, UCLA, One Stage	STT20UC


Gold UCLA Prosthetics Wide Platform Engaging		
	Size	Part #
	Restorative Kit Gold UCLA - Wide	KAUCLAOWG
	Gold UCLA Abutment, Wide (including screw)	AUCLAOWG
	Impression Transfer - Snap Wide	MPOW
	Analog – Implant Level Wide	BOW48
	Abutment Screw, UCLA, One Stage	STT20UC

Gold UCLA Prosthetics Wide Platform Non - Engaging		
	Size	Part #
	Restorative Kit Gold UCLA Non - Engaging, Wide	KAUCLAOWNG
	Gold UCLA Abutment Non-Engaging, Wide (including screw)	AUCLAOWNG
	Impression Transfer - Snap Wide	MPOW
	Analog – Implant Level Wide	BOW48
	Abutment Screw, UCLA, One Stage	STT20UC

All prices are available in the product section of our website. Volume discounts are available for some items.  
Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.  
Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann




Gold UCLA Prosthetics Narrow Platform Engaging		
	Size	Part #
	Restorative Kit Gold UCLA - Narrow	KAUCLAONG
	Gold UCLA Abutment, Narrow (including screw)	AUCLAONG
	Impression Transfer - Snap Narrow	MPON
	Analog – Implant Level Narrow	BON48
	Abutment Screw, UCLA, One Stage	STT20UC

Gold UCLA Prosthetics Narrow Platform Non - Engaging		
	Size	Part #
	Restorative Kit Gold UCLA - Narrow	KAUCLAONNG
	Gold UCLA Abutment, Narrow (including screw)	AUCLONNG
	Impression Transfer - Snap Narrow	MPON
	Analog – Implant Level Narrow	BON48
	Abutment Screw, UCLA, One Stage	STT20UC

## Angled Abutments

Only Compatible with Tissue Level implants. Angled Abutments include abutment screw. For assistance in choosing an angled abutment please see the diagram on the web site..

Gold UCLA Prosthetics Narrow Platform Engaging		
	Size	Part #
	Angled 15 degree flat oriented Regular/Narrow	AOA1565A
	Angled 15 degree corner oriented Regular/Narrow	AOA1565B

## Overdenture Abutments

Master Overdenture Restorative Kit – Chairside includes: overdenture abutment, metal housing, retention inserts, lab processing insert and block out ring.

Master Overdenture Restorative Kit – Laboratory includes: overdenture abutment, metal housing, retention inserts, lab processing insert, block out ring, analog and impression transfer.




For greater flexibility Blue Sky Bio overdenture abutments do not have an internal undercut. Instruments list above are needed for implant and abutment placement. Overdenture abutment can be placed with a .048 hex driver or a Straumann SCS driver.

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann






## Master Overdenture Restorative Kits 4.3 platform – Chairside Includes Metal Housing






	Size	Part #
	Master Overdenture Restorative Kit - Chairside, 1mm	MORK10
	Master Overdenture Restorative Kit - Chairside, 2mm	MORK20
	Master Overdenture Restorative Kit - Chairside, 3mm	MORK30

## Master Overdenture Restorative Kits 4.3 platform – Laboratory Includes Metal Housing

	Size	Part #
	Master Overdenture Restorative Kit - Laboratory, 1mm	MORK10-L
	Master Overdenture Restorative Kit - Laboratory, 2mm	MORK20-L
	Master Overdenture Restorative Kit - Laboratory, 3mm	MORK30-L

## Overdenture Products

	Size	Part #
	Complete Flat Snap Restorative Kit	NOPK
	Analog	FSA
	Replacement Inserts Strong Retention (4 Pack)	FSCS
	Replacement Inserts Medium / Light Retention (4 Pack)	FSCR
	Resin Impression transfer (2 Pack) and FSCB (4 Pack)	FSI Pack

	Size	Part #
	Processing Inserts No Retention (4 Pack)	FSCB
	Resin Impression transfer (2 Pack)	FSI
	Titanium Housing	HGFS
	.050 Hex Driver Short	HTT
	Insert Placement/Removal Instrument	TFSU

## Intraorally Preparable Abutments

Only Compatible with Tissue Level implants. Compatible with Straumann®-ITI® Lifecore® Stage 1® implants. Blue Sky Bio 1-piece angle correction abutments have increased fracture resistance and require intraoral preparation with a carbide bur and copious irrigation. Blue Sky Bio angle correction abutments are not a 2-piece design like ITI-Straumann® angled abutments.

Angle Correction Abutments will fit narrow and regular platform implants.  
Must be inserted with a .048 screwdriver head.

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

## 15° Angle Correction



Size	Part #
6mm - Preparable Abutment Regular/Narrow	AU5541

## 25° Angle Correction



Size	Part #
6mm - Preparable Abutment Regular/Narrow	AU5548

## Waxing Sleeves / Temporary Transfers


Compatible with Straumann®-ITI® Lifecore® Stage 1® implants.

Each package contains one engaging and one non-engaging sleeve, to be used for wax-ups or for temporary transfers..

### Waxing Sleeves / Temporary Transfers



Clear / Engaging




White / Non Engaging

Size	Part #
Waxing / Temporary Sleeves Narrow	WINK4100



Clear / Engaging




White / Non Engaging

Waxing / Temporary Sleeves Regular	WIRK4800
------------------------------------	----------



Clear / Engaging



White / Non Engaging

Waxing / Temporary Sleeves Wide	WIWK6500
---------------------------------	----------

## Abutment level analogs - Metal

### Abutment level analogs - Metal



Size	Part #
4mm Narrow Abutment Analog (BIO   One Stage Compatible Only)	BAN4000



5.5mm Narrow Abutment Analog (BIO   One Stage Compatible Only)	BAN5500
--	---------



Size	Part #
7mm Narrow Abutment Analog (BIO   One Stage Compatible Only)	BAN7000





4mm Regular Abutment Analog	BAR4000
-----------------------------	---------



All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

## Abutment level analogs - Metal

	Size	Part #
	5.5mm Regular Abutment Analog	BAR5500
	7mm Regular Abutment Analog	BAR7000



	Size	Part #
	4mm Wide Abutment Analog	BAW4000
	5.5mm Wide Abutment Analog	BAW5500

## Impression Transfers and Implant Level Analogs




synOcta® compatible.

Must be inserted with .048 hex drive.




### Resin Impression Transfers - Snap

	Size	Part #
	ImpressionTransfer, Narrow (BIO   One Stage Compatible Only)	MPON
	Impression Transfer, Regular	MPOR
	Impression Transfer, Wide	MPOW




### Implant level analogs

	Size	Part #
	Analog, Oct Narrow (BIO   One Stage Compatible Only)	BON41
	Analog, Oct Regular	BOR48
	Analog, Oct Wide	BOW65

### Aluminum Impression Transfers - Closed Tray

	Size	Part #
	Impression Transfer, Oct Narrow (BIO   One Stage Compatible Only)	MON5500
	Impression Transfer, Oct Regular	MOR5500
	Impression Transfer, Oct Wide	MOW5500




















### Aluminum Impression Transfers - Open Tray

	Size	Part #
	Impression Transfer, Oct Narrow (BIO   One Stage Compatible Only)	MON5500-OT
	Impression Transfer, Oct Regular	MOR5500-OT
	Impression Transfer, Oct Wide	MOW5500-OT

All prices are available in the product section of our website. Volume discounts are available for some items.  
 Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.  
 Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann



## Surgical Instruments

Surgical Instruments		
	Size	Part #
	Torque/Ratchet Wrench	TR-900-01
	Thumbwheel/ITI Wrench Adapter	ITI-Adapter
	.048 Hex Driver short	SDS048
	.048 Hex Driver long	SDL048
	J-Implant Insertion Driver short	DT0014
	J-Implant Insertion Driver long	DT0018
	Handpiece Implant Driver short	DMS
	Handpiece Implant Driver long	DML
	Drill Starter, short	DRILLTSS
	Drill Starter, long	DRILLTSL
	Drill 2.2mm, short	DRILLT22S
	Drill 2.2mm, long	DRILLT22L
	Drill 2.8mm, short	DRILLT28S
	Drill 2.8mm, long	DRILLT28L
	Drill 3.5mm, short	DRILLT35S
	Drill 3.5mm, long	DRILLT35L
	Drill 4.3mm, short	DRILLT43S
	Drill 4.3mm, long	DRILLT43L
	Drill 5.0mm, short	DRILLT50S

All prices are available in the product section of our website. Volume discounts are available for some items.  
 Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.  
 Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann



Surgical Instruments	
Size	Part #
Instrument Box	OIB

Complete Surgical Kit	CSK
-----------------------	-----



One - Stage Drill Taps	
Size	Part #
Tap 3.3mm	TAPT33
Tap 4.1mm	TAPT41
Tap 4.8mm	TAPT48
Tap 5.6mm	TAPT56



One - Stage Guide Pins	
Size	Part #
Guide Pin 2.2mm	GP0022
Guide Pin 2.8mm	GP0028
Guide Pin 3.5mm	GP0035
Guide Pin 4.3mm	GP0043



Other Surgical Instrumentation	
Size	Part #
Lindemann Bur	LINDR-S21

All prices are available in the product section of our website. Volume discounts are available for some items.  
 Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.  
 Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

## X Cube




Size	Part #
Complete Implantology System including 20:1 Handpiece (max speed 2500 RPM)	X Cube



Size	Part #
Complete Implantology and Oral Surgery System including 1:1 Straight Handpiece and 20:1 Handpiece	X Cube IOS



Size	Part #
Motor	F100aEI



Size	Part #
Back Up Implant Air Motor	SCH05-M4



Size	Part #
32:1 Handpiece External Irrigation	ACL(B)-451



Size	Part #
20:1 Handpiece External Irrigation	ACL(B)-411



Size	Part #
16:1 Handpiece External Irrigation	ACL(B)-421



Size	Part #
1:1 Straight Handpiece External Irrigation	AT-11

Size	Part #
Nozzle Clip 1:1 Straight Handpiece, Internal Irrigation	ST-Nozzle

Size	Part #
Nozzle Clip 16:1, 20:1, 32:1 Handpiece, Internal Irrigation	AG-Nozzle

Size	Part #
Y-Connector	Y-Conn



Size	Part #
Irrigation Tubing	Tube

## Sutures

All sutures are 60cm in length with reverse cutting needles and packaged 12 per box. Sutures sales are currently limited to U.S. customers.



Size	Part #
PGA Fast Resorb 4-0, FS-2 needle, 12/box	PGF40FS2



Size	Part #
PGA 4-0 violet, FS-2 needle, 12/box	PG40FS2



Size	Part #
Polypropylene (unbraided black) 5-0, P3 needle, 12/box	P50P3



Size	Part #
Nylon (unbraided black) 5-0, P3 needle, 12/box	N50P3



Size	Part #
Silk 3-0, FS-2 needle, 12/box	S30FS2

## Master Bone Expander Kit

Master Bone Expander Kit includes the Bone Expander Kit and four Lerner Ridge Splitting Chisels.

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

## Master Bone Expander Kit



Size	Part #
Master Bone Expander Kit	MBEK

## Bone Expander Kit



Size	Part #
Bone Expander Kit	BEK

## Lerner Ridge Splitting Chisels



Size	Part #
3mm Chisel with 2mm blade thickness	RSC32
3mm Chisel with .5mm blade thickness	RSC35
6mm Chisel with 2mm blade thickness	RSC62
6mm Chisel with .5mm blade thickness	RSC65

## Sinus Lift Kit

Complete Sinus Lift Kit includes: Sinus Probe, Bone Carrier, Bone Condenser, 2mm Sinus Drill, Sinus Reamers (2.4mm, 2.8mm, 3.2mm, 3.6mm), Stoppers (2mm – 11mm), Stopper Stands, and instrument box.

## Sinus Lift Kit



Size	Part #
Hammerless Atraumatic Crestal Sinus Lift Kit	Sinus Kit

## Safety Lateral Sinus Kit

## Safety Lateral Sinus Kit



Size	Part #
Safety Lateral Sinus Kit	SLSK

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann



## Terms & Conditions – Important User Information

### Blue Sky Bio, LLC - TERMS AND CONDITIONS – Important User Information

Blue Sky Bio offers high quality dental implant products at much lower prices, compared to its competitors. To help keep our overhead low and prices value-minded, we only accept returns or exchanges under *very limited circumstances*.

#### RETURNS & EXCHANGES

Blue Sky Bio will accept returns or exchanges within 30 days of original sales receipt. Products returned or exchanged between 31 and 120 days of original sales receipt are subject to a 15% restocking fee. No products can be returned or exchanged after 120 days from date of original sales receipt. Prior to return or exchange please contact us for a Return Authorization Number (RA#). All returned or exchanged products are subject to inspection. Products that have been opened, damaged or marked upon are not eligible for return or exchange. When your products arrive, inspect the packaging and the contents for damage and defects. You must notify us within 7 days from date of original sales receipt if, product packaging is damaged, illegible, or has a broken safety seal; or if product itself appears to be damaged or does not function as its intended use. In these cases, Blue Sky Bio will pay return shipping charges, and will pay shipping charges of replacement products. Blue Sky Bio does not offer a credit or reimburse for implants that do not integrate.

#### ORDER PLACEMENT

Blue Sky Bio accepts payment for purchases with Visa, MasterCard, and American Express, through its secure website or in limited circumstances, with our fax order form. Prices for products, sales tax, shipping, or any other charges will be listed on your sales receipt. Full payment is due before orders are shipped. Sales are final after shipment has left our distribution center. Blue Sky Bio reserves the right to hold orders until credit card information can be confirmed. You may be contacted by a Blue Sky Bio representative regarding your purchase.

#### SHIPPING

Orders placed within the United States are shipped via UPS. Orders placed outside the United States are shipped via Fedex. Customers are ultimately responsible to ensure that carriers will deliver at the requested time in their areas. Orders placed before 11am CST will ship the same day. Orders placed after 11am CST may not ship until the next business day. Blue Sky Bio recognizes holidays throughout the year. Watch our home page for notices regarding company closings due to holidays, as this will affect shipping delivery dates.

#### INTERNATIONAL BROKER / TAX TERMS

Blue Sky Bio utilizes Fedex Express (FE) to ship and broker all orders from the United States to your door. As broker, FE pays all taxes, tariffs, and VAT due on shipments. FE will bill you separately for these charges after you receive your order, including their broker fee. Taxes and broker fees are not part of your shipping charge collected at time of order. By placing your order you authorize Blue Sky Bio to collect any taxes not paid, as soon as reported by FE of tax delinquency, via your credit card on file. You will be responsible for any additional charges or late fees associated with not paying responsible taxes, tariffs, and VAT promptly.

#### CONTACTING BLUE SKY BIO

Contacting us is very easy and convenient. Call us toll-free 9am – 5pm CST at, **888-446-6724**, or email us anytime at, [info@blueskybio.com](mailto:info@blueskybio.com).

#### Important Product Notices For Blue Sky Bio Customers

Blue Sky Bio internal taper/oct dental implants are compatible with Straumann® implants, synOcta® system solid abutments and instrumentation. Blue Sky Bio implants have a patent-pending thread design and should not be inserted with mounts from other companies.

Straumann users require Blue Sky Bio J-implant insertion driver or handpiece driver, and the ITI Wrench Adapter. The J-implant insertion driver eliminates the need for an implant mount. All other instruments are part of the Straumann Surgical Kit.

Our implants and solid abutment system parts are compatible with Lifecore® Stage-1® implants and components.

With the exceptions listed above no other companies products should be used with our system.

Federal Law restricts these devices to sale by or on order of a licensed dentist or physician.

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

