


Compatibility Innovation Value

BIO | Trilobe



2016

Product Catalog
BIO | Trilobe

January 2016




BlueSkyBio


Table of contents


Implant Packs	2
Practice Building Kits	2,3
Cover Screws and Healing Abutments	3,4
Trilobe Conical Abutment	4
Trilobe KISS Abutments	5
Trilobe Temporary Abutments	6
Gold UCLA Abutments (including screw)	6,7
Multi Unit Abutments	7,8
Overdenture Abutments	8,9
Screw Retained Temporary Abutment	10
Implant Level Analogs	10
Implant Level Impression Transfers	10,11
Trilobe Abutment Screw	11
Surgical Instruments	11,12,13
X Cube	13,14
Safety Lateral Sinus Kit	14
Suture	14
Master Bone Expander Kit	14,15
Sinus Lift Kit	15
Terms & Condition	16


Implant Packs

Packs include Trilobe implant, 1mm cover screw and a 3mm healing abutment.
For easy insertion, implants should be placed with a J-implant insertion driver.
The patented J-implant driver prevents fracture or distortion of implants during insertion. Available only for Blue Sky Bio implants.

3.5mm Implant				
	Size	Part #		
	10mm	ITT3510		
	11.5mm	ITT3511		
	13mm	ITT3513		
	16mm	ITT3516		


4.3mm Implant				
	Size	Part #		
	8mm	ITT4308		
	10mm	ITT4310		
	11.5mm	ITT4311		
	13mm	ITT4313		
	16mm	ITT4316		


5.0mm Implant				
	Size	Part #		
	8mm	ITT5008		
	10mm	ITT5010		
	11.5mm	ITT5011		
	13mm	ITT5013		
	16mm	ITT5016		

6.0mm Implants - 5.0 Platform Switch				
	Size	Part #		
	8mm	ITT6008		
	10mm	ITT6010		
	11.5mm	ITT6011		
	13mm	ITT6013		
	16mm	ITT6016		

Practice Building Kits


Kits include Trilobe implant, 1mm cover screw, 3mm healing abutment, final abutment/implant level impression coping and implant level analog.
For easy insertion, implants should be placed with a J-implant insertion driver.
The patented J-implant driver prevents fracture or distortion of implants during insertion. Available only for Blue Sky Bio implants.

3.5mm Practice Building Kit				
	Size	Part #		
	10mm	KNPACTA3510		
	11.5mm	KNPACTA3511		
	13mm	KNPACTA3513		
	16mm	KNPACTA3516		


4.3mm Practice Building Kit				
	Size	Part #		
	8mm	KNPACTA4308		
	10mm	KNPACTA4310		
	11.5mm	KNPACTA4311		
	13mm	KNPACTA4313		
	16mm	KNPACTA4316		

All prices are available in the product section of our website. Volume discounts are available for some items.
Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.
Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

5.0mm Practice Building Kit



6.0mm Practice Building Kit - 5.0 Platform Switch



Cover Screws and Healing Abutments

Compatible with Blue Sky Bio Trilobe and NobelReplace® implants.
Cover screws and healing abutments may be placed with a .048 Hex Driver.
Blue Sky Bio cover screws and healing abutments are compatible with NobelReplace® instrumentation.

Cover Screws

	Size	Part #
	3.5mm x 1mm Cover Screw, 3.5 Platform	CLN3510
	4.3mm x 1mm Cover Screw, 4.3 Platform	CLN4310
	5.0mm x 1mm Cover Screw, 5.0 Platform	CLN5010

3mm Healing Abutments

	Size	Part #
	3.5mm x 3mm Healing Abutment	HLN3530
	4.3mm x 3mm Healing Abutment	HLN4330
	5.0mm x 3mm Healing Abutment	HLN5030
	6.0mm x 3mm Healing Abutment, 5.0 Platform Switch	HLN5030W

All prices are available in the product section of our website. Volume discounts are available for some items.
Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.
Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.


5mm Healing Abutments		
	Size	Part #
	3.5mm x 5mm Healing Abutment	HLN3550
	4.3mm x 5mm Healing Abutment	HLN4350
	5.0mm x 5mm Healing Abutment	HLN5050
	6.0mm x 5mm Healing Abutment, 5.0 Platform Switch	HLN5050W


Trilobe Conical Abutments


Only compatible with Blue Sky Bio Trilobe implants.


Abutment with included fixation screw is also impression transfer.

Conical Abutment Kits include a Conical abutment, fixation screw and analog kit contains all necessary restorative and laboratory components.

3.5 mm Abutments		
	Size	Part #
	3.5 Conical Abutment Kit	KAT3570
	3.5 Conical Abutment	ATT3570
	3.5 Implant Level Analog	BITT35

4.3 mm Abutments		
	Size	Part #
	4.3 Conical Abutment Kit	KAT4370
	4.3 Conical Abutment	ATT4370
	4.3 Implant Level Analog	BITT43

5.0 mm Abutments		
	Size	Part #
	5.0 Conical Abutment Kit	KAT5070
	5.0 Conical Abutment	ATT5070
	5.0 Implant Level Analog	BITT50

5.0 mm Abutments - Wide		
	Size	Part #
	5.0 Conical Abutment Kit	KAT5070W
	5.0 Conical Abutment	ATT5070W
	5.0 Implant Level Analog	BITT50W

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

Trilobe KISS Abutments

Only compatible with Blue Sky Bio Trilobe implants.

Complete Kits include KISS abutment, screw, metal analog, two resin analogs, impression coping, engaging and non-engaging waxing sleeves and temporary sleeve.

Impression Kits include two resin analogs, impression coping, engaging and non-engaging waxing sleeves and temporary sleeve.

3.5 Platform KISS 2mm

Size	Part #
3.5 Regular KISS Complete Kit	KASTT3502
3.5 Regular KISS Abutment 2mm collar	ASTT3502
KISS Abutment Analog, Regular	BMKISR
KISS Impression Kit Regular	KIKR

4.3 Platform KISS 2mm

Size	Part #
4.3 Regular KISS Complete Kit	KASTT4302
4.3 Regular KISS Abutment 2mm collar	ASTT4302
KISS Abutment Analog, Regular	BMKISR
KISS Impression Kit Regular	KIKR

5.0 Platform KISS 2mm

Size	Part #
5.0 Wide KISS Complete Kit	KASTT5002
5.0 Wide KISS Abutment 2mm collar	ASTT5002
KISS Abutment Analog, Wide	BMKISW
KISS Impression Kit Wide	KIKW

All prices are available in the product section of our website. Volume discounts are available for some items.


Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

Trilobe Temporary Abutments


Trilobe Temporary Abutments are only compatible with Blue Sky Bio implants.
Trilobe Temporary Abutments are not lobe engaging.

3.5 Temporary Abutment




Size	Part #
3.5 Temporary Abutment	ATTMP35

4.3 Temporary Abutment



Size	Part #
4.3 Temporary Abutment	ATTMP43

5.0 Temporary Abutment



Size	Part #
5.0 Temporary Abutment	ATTMP50

Gold UCLA Abutments (including screw)

UCLA abutments are gold alloy.

UCLA abutment includes screw.

Only compatible with Blue Sky Bio Trilobe implants. Not compatible with Nobel implants.

Gold UCLA 3.5 mm



Size	Part #
3.5 Gold UCLA Abutment	AUCLAT35G
Additional Abutment Screw	SLN18

Gold UCLA 4.3 mm



Size	Part #
4.3 Gold UCLA Abutment	AUCLAT43G
Additional Abutment Screw	SLN20





All prices are available in the product section of our website. Volume discounts are available for some items.





Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.





Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

Gold UCLA 5.0 mm		
	Size	Part #
	5.0 Gold UCLA Abutment	AUCLAT50G
	Additional Abutment Screw	SLN20









Multi Unit Abutments

3.5 Platform		
	Size	Part #
	Abutment Multi Unit Straight, 3.5 platform, Trilobe, 1.5mm cuff - Requires .080 Hex Driver	AMUSN3515
	Abutment Multi Unit Straight, 3.5 platform, Trilobe, 3.0mm cuff - Requires .080 Hex Driver	AMUSN3530
	Angled Abutment Multi Unit, 3.5 platform, 17 degree, Trilobe, includes abutment screw	AAMUN3517
	Angled Abutment Multi Unit, 3.5 platform, 30 degree, Trilobe, includes abutment screw	AAMUN3530

4.3 Platform		
	Size	Part #
	Abutment Multi Unit Straight, 4.3 platform, Trilobe, 1.5mm cuff - Requires .080 Hex Driver	AMUSN4315
	Abutment Multi Unit Straight, 4.3 platform, Trilobe, 3.0mm cuff - Requires .080 Hex Driver	AMUSN4330
	Angled Abutment Multi Unit, 4.3 platform, 17 degree, Trilobe, includes abutment screw	AAMUN4317
	Angled Abutment Multi Unit, 4.3 platform, 30 degree, Trilobe, includes abutment screw	AAMUN4330

5.0 Platform		
	Size	Part #
	Abutment Multi Unit Straight, 5.0 platform, Trilobe, 1.5mm cuff - Requires .080 Hex Driver	AMUSN5015
	Abutment Multi Unit Straight, 5.0 platform, Trilobe, 3.0mm cuff - Requires .080 Hex Driver	AMUSN5030
	Angled Abutment Multi Unit, 5.0 platform, 17 degree, Trilobe, includes abutment screw	AAMUN5017
	Angled Abutment Multi Unit, 5.0 platform, 30 degree, Trilobe, includes abutment screw	AAMUN5030

All prices are available in the product section of our website. Volume discounts are available for some items.
Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.
Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

Components					
	Size	Part #		Size	Part #
	Comfort Cap for all Multi Unit Abutments	CMU		Angled Multi Unit Abutment Delivery Aid	DSMU
	Combination Impression Transfer, Temporary Abutment and Castable Sleeve, Plastic-Requires screw SMUIHS or SMUIHL	MMUP		Screw Multi Unit Impression Transfer, long	SMUL
	Combination Impression Transfer/ Temporary Abutment - Titanium Alloy - Requires screw SMUIHS or SMUIHL	MMUT		Screw, Multi Unit, Castable Abutments, short	SMUS
	Analog, Multi Unit	BMU		.080 Hex Driver for straight Multi Unit Abutments	HFD

Overdenture Abutments

Master Overdenture Restorative Kit – Chairside includes: overdenture abutment, metal housing, retention inserts, lab processing insert and block out ring



Master Overdenture Restorative Kit – Laboratory includes: overdenture abutment, metal housing, retention inserts, lab processing insert, block out ring, analog and impression transfer

Super Snap Overdenture Restorative Kit – Chairside includes: overdenture abutment, retention inserts and block out ring

Super Snap Overdenture Restorative Kit – Laboratory includes: overdenture abutment, retention inserts, block out ring, analog and impression transfer

For greater flexibility Blue Sky Bio overdenture abutments do not have an internal undercut.

Master Overdenture Restorative Kits 3.5 platform – Chairside Includes Metal Housing

Size	Part #
 Master Overdenture Restorative Kit - Chairside, 3mm	MORKN3530
 Master Overdenture Restorative Kit - Chairside, 5mm	MORKN3550

Master Overdenture Restorative Kits 3.5 platform – Laboratory Includes Metal Housing



Size	Part #
 Master Overdenture Restorative Kit - Laboratory, 3mm	MORKN3530-L
 Master Overdenture Restorative Kit - Laboratory, 5mm	MORKN3550-L

All prices are available in the product section of our website. Volume discounts are available for some items.


Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.



Master Overdenture Restorative Kits 4.3 platform – Chairside Includes Metal Housing

	Size	Part #
	Master Overdenture Restorative Kit - Chairside, 3mm	MORKN4330
	Master Overdenture Restorative Kit - Chairside, 5mm	MORKN4350


Master Overdenture Restorative Kits 4.3 platform – Laboratory Includes Metal Housing

	Size	Part #
	Master Overdenture Restorative Kit - Laboratory, 3mm	MORKN4330-L
	Master Overdenture Restorative Kit - Laboratory, 5mm	MORKN4350-L







Master Overdenture Restorative Kits 5.0 platform – Chairside Includes Metal Housing







	Size	Part #
	Master Overdenture Restorative Kit - Chairside, 3mm	MORKN5030
	Master Overdenture Restorative Kit - Chairside, 5mm	MORKN5050

Master Overdenture Restorative Kits 5.0 platform – Laboratory Includes Metal Housing

	Size	Part #
	Master Overdenture Restorative Kit - Laboratory, 3mm	MORKN5030-L
	Master Overdenture Restorative Kit - Laboratory, 5mm	MORKN5050-L

Overdenture Products

	Size	Part #
	Complete Flat Snap Restorative Kit	NOPK
	Analog	FSA
	Replacement Inserts Strong Retention (4 Pack)	FSCS
	Replacement Inserts Medium / Light Retention (4 Pack)	FSCR
	Resin Impression transfer (2 Pack) and FSCB (4 Pack)	FSI Pack
	Ratchet Implant Driver Short	DTSBSIH

	Size	Part #
	Processing Inserts No Retention (4 Pack)	FSCB
	Resin Impression transfer (2 Pack)	FSI
	Titanium Housing	HGFS
	.050 Hex Driver Short	HTT
	.050 Hex Driver Short	SDS050
	Insert Placement/Removal Instrument	TFSU

All prices are available in the product section of our website. Volume discounts are available for some items.
Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.
Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

Screw Retained Temporary Abutment

Includes abutment screw

Base abutments can be used to bond zirconia custom abutment cylinders, angled or straight up to 15° in angulation.

3.5 Platform Abutment



Size	Part #
Ti-base for temporary abutment	AUCLAT35G
Additional Abutment Screw	SLN18

4.3 Platform Abutment



Size	Part #
Ti-base for temporary abutment	AUCLAT43G
Additional Abutment Screw	SLN20

Implant Level Analogs

Zirconia Bases



Size	Part #
3.5 Implant Level Analog	BITT35




4.3 Implant Level Analog	BITT43
--------------------------	--------




5.0 Implant Level Analog	BITT50
--------------------------	--------

Implant Level Impression Transfers



Size	Part #
3.5 Closed Tray Impression Transfer	MTT35




Size	Part #
3.5 Open Tray Impression Transfer	MTT35-OT

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.


Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.



Size	Part #
4.3 Closed Tray Impression Transfer	MTT43



Size	Part #
4.3 Open Tray Impression Transfer	MTT43-OT





Size	Part #
5.0 Closed Tray Impression Transfer	MTT50






Size	Part #
5.0 Open Tray Impression Transfer	MTT50-OT

Trilobe Abutment Screw

Only compatible with Blue Sky Bio Trilobe abutments.

Trilobe Abutment Screw		
	Size	Part #
	3.5mm Abutment Screw	SLN18
	4.3mm and 5.0mm Abutment Screw	SLN20





Surgical Instruments


Surgical Instruments		
	Size	Part #
	Adapter for NB Ratchet	NA
	Ratchet Driver	RWA
	.048 Hex Driver	SDS048

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.






Surgical Instruments		
	Size	Part #
	.048 Hex Driver	SDL048
	3.5mm J-Implant Insertion Driver short	DJT35S
	3.5mm J-Implant Insertion Driver long	DJT35L
	4.3mm J-Implant Insertion Driver short	DJT43S
	4.3mm J-Implant Insertion Driver long	DJT43L
	5.0mm J-Implant Insertion Driver short	DJT50S
	5.0mm J-Implant Insertion Driver long	DJT50L
	Starter Drill Short	DRILLNSS
	Starter Drill Long	DRILLNSL
	Drill 2.2mm Short	DRILLN22S
	Drill 2.2mm Long	DRILLN22L
	Implant Drill Short 3.5mm	DRILLN35S
	Implant Drill Long 3.5mm	DRILLN35L
	Profile Drill 3.5mm	DRILLP35
	Implant Drill Short 4.3mm	DRILLN43S
	Implant Drill Long 4.3mm	DRILLN43L
	Profile Drill 4.3mm	DRILLP43
	Implant Drill Short 5.0mm	DRILLN50S
	Implant Drill Long 5.0mm	DRILLN50L
	Profile Drill 5.0mm	DRILLP50
	Implant Drill Short 6.0mm	DRILLN60S

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

Surgical Instruments		
   	Size	Part #
	Profile Drill 6.0mm	DRILLP60
	Latch to Square Driver Converter	LS
	Trilobe Instrument Box	TIB
	Trilobe Complete Surgical Kit	TCSK

Other Surgical Instrumentation		
 	Size	Part #
	Lindemann Bur	LINDR-S21
	Paralleling Pin	GPST22

X Cube

	Size	Part #		Size	Part #
	Complete Implantology System including 20:1 Handpiece (max speed 2500 RPM)	X Cube		Complete Implantology and Oral Surgery System including 1:1 Straight Handpiece and 20:1 Handpiece	X Cube IOS
 	Size	Part #		Size	Part #
	Motor	F100aEI		1:1 Straight Handpiece External Irrigation	AT-11
	Back Up Implant Air Motor	SCH05-M4		Nozzle Clip 1:1 Straight Handpiece, Internal Irrigation	ST-Nozzle


All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

	Size	Part #
	32:1 Handpiece External Irrigation	ACL(B)-451
	20:1 Handpiece External Irrigation	ACL(B)-411
	16:1 Handpiece External Irrigation	ACL(B)-421



Size	Part #
Nozzle Clip 16:1, 20:1, 32:1 Handpiece, Internal Irrigation	AG-Nozzle
Y-Connector	Y-Conn
Irrigation Tubing	Tube

Safety Lateral Sinus Kit

Safety Lateral Sinus Kit		
	Size	Part #
	Safety Lateral Sinus Kit	SLSK

Sutures


All sutures are 60cm in length with reverse cutting needles and packaged 12 per box.
Sutures sales are currently limited to U.S. customers.


Size	Part #	Size	Part #
PGA Fast Resorb 4-0, FS-2 needle, 12/box	PGF40FS2	Nylon (unbraided black) 5-0, P3 needle, 12/box	N50P3
PGA 4-0 violet, FS-2 needle, 12/box	PG40FS2	Silk 3-0, FS-2 needle, 12/box	S30FS2
Polypropylene (unbraided black) 5-0, P3 needle, 12/box	P50P3		




Master Bone Expander Kit

Master Bone Expander Kit includes the Bone Expander Kit and four Lerner Ridge Splitting Chisels.

Master Bone Expander Kit		
	Size	Part #
	Master Bone Expander Kit	MBEK


Bone Expander Kit		
	Size	Part #
	Bone Expander Kit	BEK

All prices are available in the product section of our website. Volume discounts are available for some items.
Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.
Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

Lerner Ridge Splitting Chisels		
	Size	Part #
	3mm Chisel with 2mm blade thickness	RSC32
	3mm Chisel with .5mm blade thickness	RSC35
	6mm Chisel with 2mm blade thickness	RSC62
	6mm Chisel with .5mm blade thickness	RSC65

Sinus Lift Kit

Complete Sinus Lift Kit includes: Sinus Probe, Bone Carrier, Bone Condenser, 2mm Sinus Drill, Sinus Reamers (2.4mm, 2.8mm, 3.2mm, 3.6mm), Stoppers (2mm – 11mm), Stopper Stands, and instrument box.

Sinus Lift Kit		
	Size	Part #
	Hammerless Atraumatic Crestal Sinus Lift Kit	Sinus Kit

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

Terms & Conditions – Important User Information

Blue Sky Bio, LLC - TERMS AND CONDITIONS – Important User Information

Blue Sky Bio offers high quality dental implant products at much lower prices, compared to its competitors. To help keep our overhead low and prices value-minded, we only accept returns or exchanges under *very limited circumstances*.

RETURNS & EXCHANGES

Blue Sky Bio will accept returns or exchanges within 30 days of original sales receipt. Products returned or exchanged between 31 and 120 days of original sales receipt are subject to a 15% restocking fee. No products can be returned or exchanged after 120 days from date of original sales receipt. Prior to return or exchange please contact us for a Return Authorization Number (RA#). All returned or exchanged products are subject to inspection. Products that have been opened, damaged or marked upon are not eligible for return or exchange. When your products arrive, inspect the packaging and the contents for damage and defects. You must notify us within 7 days from date of original sales receipt if, product packaging is damaged, illegible, or has a broken safety seal; or if product itself appears to be damaged or does not function as its intended use. In these cases, Blue Sky Bio will pay return shipping charges, and will pay shipping charges of replacement products. Blue Sky Bio does not offer a credit or reimburse for implants that do not integrate.

ORDER PLACEMENT

Blue Sky Bio accepts payment for purchases with Visa, MasterCard, and American Express, through its secure website or in limited circumstances, with our fax order form. Prices for products, sales tax, shipping, or any other charges will be listed on your sales receipt. Full payment is due before orders are shipped. Sales are final after shipment has left our distribution center. Blue Sky Bio reserves the right to hold orders until credit card information can be confirmed. You may be contacted by a Blue Sky Bio representative regarding your purchase.

SHIPPING

Orders placed within the United States are shipped via UPS. Orders placed outside the United States are shipped via Fedex. Customers are ultimately responsible to ensure that carriers will deliver at the requested time in their areas. Orders placed before 11am CST will ship the same day. Orders placed after 11am CST may not ship until the next business day. Blue Sky Bio recognizes holidays throughout the year. Watch our home page for notices regarding company closings due to holidays, as this will affect shipping delivery dates.

INTERNATIONAL BROKER / TAX TERMS

Blue Sky Bio utilizes Fedex Express (FE) to ship and broker all orders from the United States to your door. As broker, FE pays all taxes, tariffs, and VAT due on shipments. FE will bill you separately for these charges after you receive your order, including their broker fee. Taxes and broker fees are not part of your shipping charge collected at time of order. By placing your order you authorize Blue Sky Bio to collect any taxes not paid, as soon as reported by FE of tax delinquency, via your credit card on file. You will be responsible for any additional charges or late fees associated with not paying responsible taxes, tariffs, and VAT promptly.

CONTACTING BLUE SKY BIO

Contacting us is very easy and convenient. Call us toll-free 9am – 5pm CST at, **888-446-6724**, or email us anytime at, info@blueskybio.com.

IMPORTANT PRODUCT NOTICES FOR BLUE SKY BIO CUSTOMERS

Blue Sky Bio Two Stage Trilobe implants are compatible with NobelReplace® prosthetics and instrumentation. Blue Sky Bio implants have a patented driver system that prevents fracture during insertion. The driver system is compatible only with Blue Sky Bio implants and cannot be used with any other systems. Only use Blue Sky Bio smooth tapered drills.

Blue Sky Bio's advanced Trilobe conical connection abutments are to be used only with Blue Sky Bio implants.

With the exceptions listed above no other companies products should be used with our system.

Federal Law restricts these devices to sale by or on order of a licensed dentist or physician.

Blue Sky Bio is not affiliated with Nobel Biocare. NobelReplace® is a registered trademark of Nobel Biocare.

Blue Sky Bio tri-lobe implants and abutments are sold subject to a conditional license under U.S. Patent No. 6,733,291 and U.S. Patent No. 5,795,160. Purchasers of Blue Sky Bio tri-lobe implants and abutments agree that the conditional license does not extend to the use of such products in conjunction with products not manufactured by, or on behalf of, Blue Sky Bio or Nobel Biocare.

